



US006058820A

# United States Patent [19] Rinner

[11] **Patent Number:** **6,058,820**  
[45] **Date of Patent:** **\*May 9, 2000**

- [54] **ROD CUTTER**
- [75] Inventor: **James A. Rinner**, Racine, Wis.
- [73] Assignee: **Beere Precision Medical Instruments, Inc.**, Racine, Wis.
- [\*] Notice: This patent issued on a continued prosecution application filed under 37 CFR 1.53(d), and is subject to the twenty year patent term provisions of 35 U.S.C. 154(a)(2).

812,452	2/1906	Richards	83/200
1,265,345	5/1918	Rock	83/200 X
1,874,616	8/1932	Porter	30/193 X
2,527,735	10/1950	Johnson	30/93 X
2,649,913	8/1953	Linder	83/200
2,910,900	11/1959	Klein	30/193 X
5,261,303	11/1993	Strippgen	83/199

### FOREIGN PATENT DOCUMENTS

372	2/1869	United Kingdom	83/199
4029	2/1907	United Kingdom	83/199

- [21] Appl. No.: **09/069,697**
- [22] Filed: **Apr. 29, 1998**
- [51] Int. Cl.<sup>7</sup> ..... **B26D 3/16**
- [52] U.S. Cl. .... **83/200; 83/633; 30/93; 30/253**
- [58] Field of Search ..... 30/253, 191, 193, 30/175, 176, 186, 93, 94, 95; 83/633, 634, 199, 200

*Primary Examiner*—Rinaldi I. Rada  
*Assistant Examiner*—Charles Goodman  
*Attorney, Agent, or Firm*—Arthur J. Hansmann

### [57] ABSTRACT

A rod cutter having two force-applying members pivotally attached together. A cutter is pivotally mounted on each of the members in a mechanical advantage arrangement. The cutters both have extending ends which are pivotally connected together and have rod cutting edges thereon, all for force-advantage cutting.

- [56] **References Cited**  
U.S. PATENT DOCUMENTS  
690,083 12/1901 Stolpe ..... 83/199

**21 Claims, 3 Drawing Sheets**

