



US006739068B1

(12) **United States Patent**
Rinner

(10) **Patent No.:** **US 6,739,068 B1**
(45) **Date of Patent:** **May 25, 2004**

- (54) **PLIERS WITH JAW SPACING AND LOAD MEASURING READINGS**
- (75) Inventor: **James A. Rinner, Racine, WI (US)**
- (73) Assignee: **Pilling Weck Incorporated, Horsham, PA (US)**
- (*) Notice: Subject to any disclaimer, the term of this patent is extended or adjusted under 35 U.S.C. 154(b) by 0 days.
- (21) Appl. No.: **10/336,677**
- (22) Filed: **Jan. 6, 2003**
- (51) **Int. Cl.⁷** **G01B 5/02; A61B 17/66**
- (52) **U.S. Cl.** **33/783; 33/806**
- (58) **Field of Search** 33/783, 784, 792, 33/794, 797, 798, 800, 801, 806, 807, 813, 815, 512, 542; 7/125; 606/57

4,127,112 A	11/1978	Sherlock et al.	128/25
4,184,259 A	1/1980	Sosnay	433/4
4,226,025 A	10/1980	Wheeler	33/148
4,312,363 A	1/1982	Rothfuss et al.	128/774
4,432,376 A *	2/1984	Huszar	33/784
4,633,587 A	1/1987	Harrison	30/261
4,843,721 A	7/1989	Hoge	33/802
4,898,161 A *	2/1990	Grundeis	606/105
5,070,623 A	12/1991	Barnes	33/807
5,116,340 A	5/1992	Songer et al.	606/103
5,156,161 A	10/1992	Lollar	128/774
5,183,055 A	2/1993	Seager	
5,297,538 A *	3/1994	Daniel	600/206
5,336,228 A	8/1994	Cholhan	606/119
5,381,799 A	1/1995	Hamilton et al.	128/777
5,433,015 A *	7/1995	Mazenet	33/815
5,597,305 A	1/1997	Ray	
5,658,295 A	8/1997	Krementsov	606/119
6,048,322 A	4/2000	Kushida	600/587
6,102,909 A *	8/2000	Chen et al.	606/45
6,263,770 B1	7/2001	Gomas et al.	81/427
6,551,316 B1 *	4/2003	Rinner et al.	606/57
6,553,685 B2 *	4/2003	Nishina et al.	33/815

(56) **References Cited**

U.S. PATENT DOCUMENTS

466,986 A	1/1892	Van Roden	
1,528,273 A	3/1925	Shwed	
1,626,540 A	4/1927	Kimura	
2,112,873 A	4/1938	Wright	81/18
2,125,945 A	8/1938	Montgomery	
2,396,619 A	3/1946	Strayer	
2,591,360 A	4/1952	Keenan	81/47
2,888,825 A	6/1959	Kraft	
3,140,546 A	7/1964	Bartlett	33/148
3,140,715 A	7/1964	Whitton	
3,391,573 A	7/1968	Hiller	
3,559,515 A	2/1971	Kyser	
3,740,779 A	6/1973	Rubricvis	7/14.1 R
4,050,464 A	9/1977	Hall	

* cited by examiner

Primary Examiner—Christopher W. Fulton
(74) *Attorney, Agent, or Firm*—Arthur J. Hansmann

(57) **ABSTRACT**

A pliers with two handles and two jaw portions, all pivoted together. Two readable combinations, each having a scale and a pointer, are on the handles to respectively reveal the amount of pivot and the force applied to a workpiece. A parallel linkage assembly can be attached to the jaw portions for contacting the workpiece, and, because of the parallelism, the location of the contact with the workpiece is not critical to the readings.

23 Claims, 4 Drawing Sheets

